

**HAORUN® MEDICAL**



HAORUN MEDICAL DRESSING CO.,LTD  
www.haorunmed.com sales@haorunmed.com  
Mobile/Whatsapp/Wechat: +86 13566542343



**HAORUN MEDICAL DRESSING CO.,LTD**



# COMPANY PROFILE

HAORUN Medical Dressing Co., Ltd (Group) is both a producer and trader of Medical wound dressings, Operating room supplies, Respiratory and Urinary products, Electronic, Monitoring products, Household hygiene products as well as cosmetics. After years of development, HAORUN now has its own workshops for manufacturing and packaging as well as efficient teams for procurement, quality control and sales. The quality management system of HAORUN has successfully passed ISO 13485:2016 (TUV). The products of HAORUN has reached the various local standards of China, U.S., E.U., etc, and it has obtained CE and FSC and other certificates.

HAORUN has won the trust of numerous customers from over 90 countries all over South Americas, Europe, the Middle East, Africa, South East Asia as well as Mainland China. To better serve its customers, HAORUN may offer OEM service.

HAORUN attaches great importance to product quality and customers' needs. We will always provide the best service for its customers with a high sense of public quality and social responsibility. HAORUN wishes to cooperate with more customers and to make the world a better place together.



America	Europe
Bolivia	Belgium
Chile	Bulgaria
Colombia	France
Costa Rica	Greece
Dominica	Lithuania
Haiti	Italy
Honduras	Poland
Jamaica	Portugal
Peru	Romania
Uruguay	Serbia
Venezuela	Spain
	Ukraine

Africa	
Algeria	Rwanda
Botswana	Senegal
Congo	Sudan
Egypt	South Africa
Eritrea	Tanzania
Ethiopia	Uganda
Ghana	Zambia
Kenya	Zimbabwe
Libya	
Mautitius	
Moldova	
Mozambique	

Asian		
Byanmar	Bahrain	Alghanistan
Bengal	Qatar	India
Indonesia	Lebanon	Kazakhstan
Laos	Palestine	Maldives
Malaysia	Saudi Arabia	Pakistan
Philippines	Turkey	Sri Lanka
Singapore	Yemen	Uzbekistan
Thailand	Iraq	
Vietnam	Jordan	
	United Arab Emirates	

Contact us :

HAORUN MEDICAL DRESSING CO.,LTD  
 www.haorunmed.com sales@haorunmed.com  
 Mobile/Whatsapp/Wechat: +86 13566542343

CONTENTS

01-05 Gauze Series

06-10 Non-woven Series

11-14 Cotton Series

15 Notes



**HAORUN<sup>®</sup> MEDICAL**

**Gauze Series**

- Absorbent Gauze Roll
- Pillow Gauze
- Gauze Swabs
- Gauze Sheet
- Cutting Gauze
- Laparotomy Sponges
- Gauze Bandages
- Lint Gauze Roll
- Triangular Bandages
- Fluffy Bandages
- Gauze Ball
- Vaseline Gauze



## Absorbent Gauze Roll

100% pure cotton  
BP 11/13/17/20 threads  
white/green/blue color  
with/without X-ray  
1/2/4-ply

\*can be produced as requirements



## Pillow Gauze

100% pure cotton  
BP 11/13/17/20 threads  
white/green/blue color  
with/without X-ray

\*can be produced as requirements



## Gauze Swabs

100% pure cotton BP 11/13/17/20 threads  
white/green/blue color  
with/without X-ray folded/unfolded edges  
5X5/7.5X7.5/10X10/10X15/10X20cm etc.  
8/12/16/24/32-ply  
non-sterile/sterile

\*can be produced as requirements



## Gauze Sheet

100% pure cotton BP 11/13/17/20 threads  
white/green/blue color  
36"X100yards 44"X22yards etc.

\*can be produced as requirements



## Cutting Gauze

100% pure cotton BP 11/13/17/20 threads  
white/green/blue color  
with/without X-ray  
10X10/10X20/20X20/30X30/40X40cm etc.  
non-sterile/sterile


\*can be produced as requirements



## Laparotomy Sponges

100% pure cotton BP 11/13/17/20 threads  
white/green/blue color  
with/without X-ray or chip  
with/without cotton loop  
10X60/20X20/30X30/45X45/50X50cm etc.  
4/6/8/12-ply  
non-sterile/sterile unwashed/pre-washed

\*can be produced as requirements



## W.O.W Gauze Bandages

100% pure cotton BP 11/13/17/20 threads  
cutting/woven edges  
width:5/7.5/10/12.5/15/20cm etc.  
length:4.5/5/10m etc.  
12rolls per dozen/individual bag

\*can be produced as requirements



## Lint Gauze Roll

100% pure cotton  
454-500g etc.

\*can be produced as requirements



## Triangular Bandages

100% pure cotton  
unbleached  
90X90X131/90X90X127cm etc.  
safety pin

\*can be produced as requirements



## Fluffy Bandages

100% pure cotton BP 13/17/20 threads  
width:5/7.5/10cm etc.  
length:10m etc.  
non-sterile/sterile

\*can be produced as requirements



## Gauze Ball

100% pure cotton BP 17/20 threads  
with/without X-ray  
8X8/10X10/15X15/20X20/30X30/40X40cm  
non-sterile/sterile

\*can be produced as requirements



## Vaseline Gauze

100% pure cotton  
21S/24X24 200±20g/m<sup>2</sup>  
10X10/20X20/30X30/40X40cm etc.  
individual aluminum foil bag/plastic box  
gama-sterile

\*can be produced as requirements



# HAORUN<sup>®</sup> MEDICAL

## Non-woven Series

- Non-woven Swabs
- Drain Sponges
- Non-woven Ball
- Abdominal Pad (ABD)
- Underpads
- Adult Diapers
- Isolation Gown
- Surgical Gown
- Disposable Drawsheet
- Surgical Drapes
- Non-woven Cap
- Non-woven Shoes Cover



### Non-woven Swabs

70% viscose + 30% polyester  
 30/35/45/50gsm  
 with/without X-ray  
 5X5/7.5X7.5/10X10cm etc.  
 4/6/8-ply  
 non-sterile/sterile

\*can be produced as requirements



### Drain Sponges

70% viscose + 30% polyester  
 30/35/45/50gsm  
 Y-cut/I-cut  
 5X5/7.5X7.5/10X10cm etc.  
 4/6/8-ply

\*can be produced as requirements



### Non-woven Ball

100% pure cotton  
 30/35/45/50gsm  
 with/without X-ray  
 10X10/20X20/25X25/30X30/30X40cm etc.  
 non-sterile/sterile


\*can be produced as requirements



**Abdominal Pad (ABD)**

non-woven cover outside  
tissue and cellulose inside  
all edges are sealed to prevent linting  
5"X9" 7.5"X8" 8"X10" etc.  
non-sterile/sterile

\*can be produced as requirements



**Underpads**

top layer: non-woven fabric and PE film  
middle layer: consists of imported fluff and SAP  
60X40/60X45/60X60/75X45/75X60/90X60  
80X150cm etc.

\*can be produced as requirements



**Adult Diapers**

non-woven fabric surface  
tissue paper+SAP+fluff pulp+blue ADL  
PP/magic tape  
elastic waistband  
M/L/XL/XXL

\*can be produced as requirements



**Isolation Gown**

SMS/PP+PE/CPE/SBPP/SFS  
18-50gsm  
blue/white/green/yellow color  
short/long sleeve  
S/M/L/XL

\*can be produced as requirements



**Surgical Gown**

SMS/PP+PE/CPE/SBPP/SFS  
40-50gsm  
blue/white/green/yellow color  
short/long sleeve  
M/L/XL/XXL

\*can be produced as requirements



**Disposable Drawsheet**

paper and PE film  
80X210cm 100X200cm etc.

\*can be produced as requirements



**Surgical Drapes**

non-woven/SMS  
40-60gsm  
blue/green  
40X60/50X70/60X80cm

\*can be produced as requirements



**Non-woven Cap**

PP non-woven  
clip/bouffant/doctor cap  
single/double elastic  
white/blue/green/pink/yellow color

\*can be produced as requirements



**Non-woven Shoes Cover**

PP non-woven  
single/double elastic  
blue/green/pink color  
15X40/15X36cm etc.

\*can be produced as requirements



**HAORUN<sup>®</sup> MEDICAL**

**Cotton Series**

- Absorbent Cotton Wool
- Absorbent Cotton Wool
- Zig-zag Cotton Wool
- Cotton Ball
- Cotton Pad
- Cotton Tipped Applicator
- Absorbent Dressing Pad
- Cosmetic Cotton
- Dental Cotton Roll






## Absorbent Cotton Wool

100% pure cotton  
50/100/200/250/400/454/500/1000kg

can be produced as requirements



## Absorbent Cotton Wool

100% pure cotton  
50/100/200/250/400/454/500/1000kg  
interleaved paper  
cutting sides

\*can be produced as requirements



## Zig-Zag Cotton Wool

100% pure cotton  
50/100/200/250/400/454/500/1000kg

\*can be produced as requirements



## Cotton Ball

100% pure cotton  
0.2/0.3/0.4/0.5/0.8/1.0/2.0g  
non-sterile/sterile

\*can be produced as requirements



## Cotton Pad

100% pure cotton  
different size/gram weight  
4X4/5X5/7.5X7.5/10X10cm  
non-sterile/sterile

\*can be produced as requirements



## Cotton Tipped Applicator

100% pure cotton  
wooden/plastic stick  
single/double tip  
different size  
non-sterile/sterile

\*can be produced as requirements



**Absorbent Dressing Pad**

100% cotton  
gauze cover  
cutting/sealed/sewed edges  
non-sterile/sterile  
7X10/10X10/10X15/20X20/25X45cm etc.

\*can be produced as requirements



**Cosmetic cotton**

100% pure cotton  
square/round type  
0.4/0.55/0.6g

\*can be produced as requirements



**Dental Cotton Roll**

size #1 0.8X3.8cm  
size #2 1.0X3.8cm  
size #2 1.2X3.8cm  
size #2 1.5X3.8cm  
non-sterile/sterile

\*can be produced as requirements

Series of horizontal dashed lines for taking notes.